Conditional Diffusion for EEG Data Synthesis


Semester Project for Advanced Machine Learning 2024

27.05.2024


Group 18

Introduction


- EEG: the brain activity of different areas of the cerebral outer cortices
- Dataset (TUH EEG) with labels for different kinds of seizure activity
- Generally, seizures are slow waves between 2 10 Hz, hyper-synchronous in time
- Idea: Generate synthetic EEG with seizure and non-seizure activity




Healthy EEG



Seizure EEG




Overview




- 1. Train a diffusion model using real EEGs with two targets (seizure, non-seizure)
- 2. Use diffusion model to sample synthetic EEGs (conditioned by target)
- 3. Use synthetic EEGs to train a (simple) classifier model
- 4. Evaluate classifier model on real EEGs

Model U-Net with Attention



Results



Validation


(w/ perfect class balance)

Accuracy: 0.80 TPR: 0.89 TNR: 0.71


F1 score: 0.81 AUROC: 0.85

Accuracy: 0.5 TPR: 1.0 TNR: 0.0


F1 score: 0.66 AUROC: 0.70




a) Train: Real - Test: Real



c) Train: Real - Test: Synth



b) Train: Synth - Test: Synth



c) Train: Synth - Test: Real

Accuracy: 0.57 TPR: 0.30 TNR: 0.84 F1 score: 0.41 AUROC: 0.77

Accuracy: 1.0

F1 score: 1.0

AUROC: 1.0


TPR: 1.0

TNR: 1.0

Thank you for your attention

Questions?

Appendix A




a) Real samples

b) Synthetic samples

Frequency features

Appendix B





a) Real samples

b) Synthetic samples


Frequency features


Appendix C: Channels and classes





Appendix D: Original signal and spectrum






Appendix E: Choice of classifier

- Try out different classifier models
- Train on real data, test on real data
- Features:
 - o FFT
 - Discrete wavelet transform
 - Mean, Standard deviations, Percentiles
 - Welch's PSD
- Models:
 - Gradient boosting
 - o SVC
 - Multilayer perceptron

Model chosen: FFT + SVC

Appendix F



Appendix G: Dataset

- https://isip.piconepress.com/projects/tuh_eeg/html/downloads.shtml
- Access on request
- 26 846 clinical EEG recordings collected at Temple University Hospital (TUH) from 2002 - 2017
- EEGs stored as EDF files (EDF = European Data Format)
- Annotations stored as CSV (channels, start, stop, class [e.g. fnsz])